


## Low Temperature Constant Temperature Bath-Vertical Water Tank

- Vertical Constant Temperature Low Temperature Tank PID Intelligent temperature control, temperature control timing, heating and cooling and other functions.
- It can be operated in constant temperature refrigeration, internal circulation and external circulation modes. The external circulation can be operated with other instruments.
- The temperature control range can reach  $-30\sim 100^{\circ}\text{C}$ , and other ranges can be customized. LCD The digital display screen can display the temperature and timing time in real time.
- There is a drain valve on the side of the instrument to facilitate the discharge of liquid.
- The instrument adopts iron spray-painted shell, stainless steel liner and table. The stainless steel liner is acid and alkali resistant, corrosion resistant and easy to clean.

*\* Use clean water for constant temperature medium above  $5^{\circ}\text{C}$ .*

*\* Use 75% or more alcohol below  $5^{\circ}\text{C}$*

model	DC-0006	DC-0506	DC-1006	DC-2006	DC-3006	Image
Temperature range	0~100°C	-5~100°C	-10~100°C	-20~100°C	-30~100°C	
Accuracy	$\pm 1^{\circ}\text{C}$					
Volatility	$\pm 0.2^{\circ}\text{C}$					
Resolution	0.7					
Cycle mode	Internal/External Loop					
Timing duration	9999 minutes					
Pump flow	6.5L/min					
Lift	1.5m					
Tank capacity	6L	10L	20L	30L	50L	
Slot opening	185*150MM	185*150MM	260*220	315*310MM	315*310MM	
Groove Depth	150MM	150MM	150MM	200MM	200MM	
Drainage function	Yes	Yes	Yes	Yes	Yes	
Dimensions	402*263*392	405*263*443	483*355*438	583*405*438	/	